

A NEW RANGE OF HIGH PERFORMANCE, HEAVY DUTY WHEELCHAIRS
ENGINEERED TO MEET ALL YOUR REQUIREMENTS.



Moly - lightweight and versatile



Phönix - versatile and strong



Goliath - strong and safe

MANUFACTURED BY

SUMED[®]

Moly - Technical specification (dimensions in cm)

Seat Width	Seat Height	Seat Depth	Backrest Height	Total Width	Total W. folded	Total Height ¹	Total Length ²	Armrest Height	Transport Weight ³	Total Weight	Max. Load
36	42.5-52.5	40-44	36-46	56/58	30	95	103	22/27	11	20 _{approx.}	125
38	42.5-52.5	40-44	36-46	58/60	30	95	103	22/27	11	20 _{approx.}	125
40	42.5-52.5	40-44	36-46	60/62	30	95	103	22/27	11	21 _{approx.}	125
43	42.5-52.5	40-44	36-46	63/65	30	95	103	22/27	11	21 _{approx.}	125
45	42.5-52.5	40-44	36-46	65/67	30	95	103	22/27	11	22 _{approx.}	125
48	42.5-52.5	40-44	36-46	68/70	30	95	103	22/27	11	22 _{approx.}	125
51	42.5-52.5	40-44	36-46	71/73	30	95	103	22/27	11	23 _{approx.}	125
53	42.5-52.5	40-44	36-46	73/75	30	95	103	22/27	11	26 _{approx.}	150
56	42.5-52.5	40-44	36-46	76/78	30	95	103	22/27	11	26 _{approx.}	150

¹Total width without/with drum brake ²Total Height with Standard Seat Height 50 cm and Standard Backrest Height 42 cm ³Total length without wheelbase elongation

⁴Transport Weight without rear wheels, armrests and legrests

Phönix - Technical specification (dimensions in cm)

Seat Width	Seat Height	Seat Depth	Backrest Height	Total Width	Total W. folded	Total Height ¹	Total Length	Armrest Height	Transport Weight ²	Total Weight	Max. Load
50	42.5-52.5	42-50	36-46	70	29	92	108	14/24	14	25 _{approx.}	185
53	42.5-52.5	42-50	36-46	73	29	92	108	14/24	14	25 _{approx.}	185
56	42.5-52.5	42-50	36-46	76	29	92	108	14/24	14	25 _{approx.}	185
59	42.5-52.5	42-50	36-46	79	29	92	108	14/24	14	25 _{approx.}	185
62	42.5-52.5	42-50	36-46	82	29	92	108	14/24	14	25 _{approx.}	185
65	42.5-52.5	42-50	36-46	85	29	92	108	14/24	14	25 _{approx.}	185

¹Total Height with Standard Seat Height 50 cm and Standard Backrest Height 42 cm ²Transport Weight without rear wheels, armrests and leg rests

Goliath - Technical specification (dimensions in cm)

Seat Width	Seat Height	Seat Depth	Backrest Height	Total Width	Total W. folded	Total Height ¹	Total Length	Armrest Height	Transport Weight ²	Total Weight	Max. Load
55	45-54	45-54	42or46	75	38	94	104	15/25	21	41 _{approx.}	300
60	45-54	45-54	42or46	80	38	94	104	15/25	21	41 _{approx.}	300
65	45-54	45-54	42or46	85	38	94	104	15/25	21	41 _{approx.}	300
70	45-54	45-54	42or46	90	38	94	104	15/25	21	41 _{approx.}	300
75	45-54	45-54	42or46	95	38	94	104	15/25	21	41 _{approx.}	300
80	45-54	45-54	42or46	100	38	94	104	15/25	21	41 _{approx.}	300
85	45-54	45-54	42or46	105	38	94	104	15/25	21	41 _{approx.}	300

¹Total Height with Standard Seat Height 51 cm and Standard Backrest Height 42 cm ²Transport Weight without rear wheels, armrests and leg rests



Sumed, Uniroll and Rollo are trademarks of Sumed.

MANUFACTURED BY

SUMED[®]

Sumed International (UK) Ltd.
Integrity House, Units 1-2 Graphite Way, Hadfield, Glossop, Derbyshire. SK13 1QH
Tel: +44 (0)1457 890980 Fax: +44 (0)1457 890990
Email: sales@sumed.co.uk www.sumed.co.uk