

Community Mattresses from **SUMED**[®]

Solutions for Medium and High Risk Clients

Viscotech[®] Community range has been specifically designed to deliver excellent pressure relief at a very competitive price and could offer you real clinical advantages and value for money - with prices from as low as **£110 per mattress**.

Dimensions 197 x 88 x 15 cm

All Viscotech[®] Community Mattresses are manufactured with a high quality polyurethane cover incorporating welded seams and a full depth overlap to prevent fluid ingress. The white interior helps to highlight any fluid penetration.

Two variants are available for Viscotech[®] Community Mattresses. The first is a Medium Risk Mattress with a maximum user weight limit of 28 Stones (178kg) double castellated at a weight of 10.3kg.

There is also a High Risk Mattress available with a maximum user weight limit of 40 Stones (254kg) non-turn and single castellated. The Mattress including the cover are at a weight of 16kg.

The Viscotech[®] Community Mattress is made from a high quality durable 6" (15cm) castellated foam. The specially contoured castellations are smaller under the head and heel areas and larger under the sacral area for maximum patient immersion and comfort.

Call us now to evaluate Viscotech[®] Community Mattress for yourself.

Please contact our dedicated team for more information:


Sumed International (UK) Ltd
Integrity House
Units 1-2 Graphite Way
Hadfield, Glossop
Derbyshire SK13 1QH
Phone: 01457 890980
Fax: 01457 890990
E-mail: sales@sumed.co.uk
www.sumed.co.uk



Viscotech[®] Community - High Risk



Viscotech[®] Community - Medium Risk



XSensor[®] Pressure Mapping of Viscotech[®] Community Medium Risk Mattress

Suggested Companion Products



Viscotech[®] Cushions are available up to Very High Risk and 28 Stones (178kg) weight limit